Epicentre QuickExtract Family 相關常見問題

1. Epicentre QuickExtract Family 有哪些？適用於哪些樣品？

產品名稱	產品編號	標準操作時間	後續應用	適用樣品
QuickExtract™ DNA Extraction Solution	QE0905T (5ml) QE09050 (50ml)	65°C 6min →98°C 2min (difficult lysis procedure: 65°C 30min → 98°C 15min)	Genotyping, genetic studies, identity testing, viral/microbial screening	Hair follicles, feathers (quill-end cells), tissue-culture cells, buccal cells, zebrafish (organs, scales), mouse tail snips
QuickExtract™ RNA Extraction Kit	QER09015 (5ml) QER090150 (50ml)	1,000 x g 4°C 5min remove supernatant → +100ul solution vortex 1min	RT-PCR	Cultured cells: human, mouse, rat, <i>E. coli</i> , <i>S. aureus</i>
QuickExtract™ Plant DNA Extraction Solution	QEP80705 (5ml) QEP70750 (50ml)	65°C 6min →98°C 2min	PCR, e.g., GMO testing	<i>Arabidopsis</i> , barley, maize, emmer, pepper, rice, spelt, spinach, soybeans, wheat
QuickExtract™ Seed DNA Extraction Solution	QES08095T (5ml) QES080950 (50ml)	65°C 6min →98°C 2min	PCR, e.g., GMO testing	Apple, cotton, sunflower, tomato, barley, maize, oats, rice, rye, wheat
QuickExtract™ Bacterial DNA Extraction Kit	QEB0905T (5ml) QEB09050 (50ml)	RT 15min →80°C 2min	PCR, restriction digests, PFGE, optical mapping	Gram-positive: <i>Bacillus subtilis</i> , <i>Bifidobacterium spp.</i> , <i>Brevibacterium linens</i> , <i>Clostridium perfringens</i> , <i>Lactobacillus plantarum</i> , <i>Listeria monocytogenes</i> , <i>Staphylococcus equorum</i> , <i>Streptococcus agalactiae</i> , <i>Streptococcus pyrogenes</i> , <i>Streptococcus thermophilus</i> Gram-negative: <i>E. coli</i> , <i>Salmonella typhimurium</i> , <i>Vibrio gazogenes</i>
QuickExtract™ FFPE DNA Extraction Kit	QEF81805 (5ml) QEF81050 (50ml)	56°C 60min →98°C 2min	Microsatellite, single-nucleotide polymorphisms (SNP), tumor heterogeneity studies, copy number variations (CNV), methylation analysis, short tandem repeats (STR)	Formalin-fixed, paraffin-embedded human tissue samples
QuickExtract™ FFPE RNA Extraction Kit	QFR82805 (5ml) QFR82050 (50ml)	56°C 30min →98°C 2min	RT-PCR	Formalin-fixed, paraffin-embedded human tissue samples

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4. QuickExtract DNA Extraction solution 何時該使用 difficult lysis procedure?

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