Composition for Maintaining Organ and Cell Viability

- ⇒ Organ Preservation and Transportation
- ⇒ Cell Culture, FBS Replacement
- ⇒ Tissue Preservation & Transportation

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- ◆ Animal/Human Component Free
- ◆ Low Potassium
- ◆ Dextran Based
- ◆ Liposome Complex
- ◆ Nutrients
- ◆ GMP Validated Manufacturing
- ◆ Batch Consistency Lot-to-Lot
- ◆ US Patent No 7,220,538

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Composition of Lifor®: Physiological, oxygen enriched solution, inorganic salts, amino acids, vitamins, adenosine, cholesterol, glucose, dextrane 70, growth factors (EGF, VEGF, HGF).

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
Figure 1: 保存在 Lifor® 中不同的保存時間及溫度的癌組織抽取 RNA,皆沒有 RNA 降解的現象 : **A.** PC-3 cells for up to 72h at RT, **B.** Lung Cancer tissue 5 days at RT and 4°C, respectively and **C.** PC-3 cells for up to 3 months at -196°C. RNA did not show any degradation.

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10h Preservation of Guinea Pig Isolated Hearts Perfused at Low Flow with Air-saturated Lifor® solution at 26°C: comparison to Viaspan®, David F. Stoe, Amadou K.S. Camara, James S. Heisner, Mohammed Aldakkak, David R. Harder, American Physiological Society, 2007, vol. 62

Preservation of Hearts with Air-saturated Lifor Solution at Room Temperature for 20 Hours, D.F. Stowe, M. Aldakkak, J.S. Heisner, AKS Camara, D.R. Harde, *The Journal of Heart and Lung Transplant*, September 2008

Preservation of Hearts with Air-saturated Lifor Solution at Room Temperature for 20 Hours, D.F. Stowe *et al*, *The Journal of Heart and Lung Transplant*, September 2008

A Novel Method for Generating Xeno-Free Human Feeder Cells for Human Embryonic Stem Cell Culture, G. Meng, *et al.*, *Stem Cells and Development*, October 2007

Cellular Incorporation Into Electrospun Nanofibers Retained Viability, Proliferation, and Function in Fibroblasts, John A. van Aalst, *et al.*, SOUTHEASTERN SOCIETY OF PLASTIC AND RECONSTRUCTIVE SURGEONS, 2008

Establishing continuous single cell cultures from primary tumor excisions using as new generation of tumor tissue transport- and culture solution, R.A. Hilger *et al.*,

Nonfrozen Transport Medium Preserves and Restores Skeletal Muscle Enzymatic Activity and Morphology, Iren Horkayne-Szakaly, *et al.*, *The Journal of Histotechnology*, Vol. 32 No. 2, June 2009.

A comparative study of freezing single cells and spheroids: Towards a new model system for optimizing freezing protocols for cryobanking of human tumours
F. Ehrhart, *et al.*, *Cryobiology*, 2008.

Establishment of a transport system for mouse epididymal sperm at refrigerated temperatures. Takeo T, *et al.*, *Cryobiology*. 2012 Jun 18.