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
TUNEL-based DNA Fragmentation Assay Kits

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


2. Dilute wash buffer 1:10 with diH₂O.

3. Further dilute FLICA[®] with 200 μ L PBS.

4. Add 5-10 μ L FLICA[®] to each sample (~300 μ L aliquot of cultured cells).

5. Incubate 15-60 minutes.

6. Remove media.

7. Wash cells: add wash buffer and spin cells (twice); or add fresh media and incubate 30 minutes (twice).


8. If desired, label with additional stains, such as Hoechst, Propidium Iodide, 7-AAD, or an antibody.

9. If desired, fix or embed cells.


10. Analyze with a fluorescence microscope, plate reader, or flow cytometer.

Apoptosis Inducer	Cat. No. and Size	Note
Camptothecin	6208 (25mg), 6209 (100mg), 6210 (250mg)	It is a prototypic DNA topoisomerase I inhibitor
Staurosporine	6212 (1mg)	It induces DNA fragmentation and apoptosis
Granzyme B	6211 (5000U)	It is a component of the cytolytic granules of cytotoxic T lymphocytes and natural killer cells
TNF-alpha	6213 (100 μ g)	It is purified from recombinant E. Coli as a single, non-glycosylated peptide containing 158 amino acids with a MW of 17.5 kDa.
TNF-beta	6214 (20 μ g)	It is purified from recombinant E. Coli as a single, non-glycosylated peptide containing 171 amino acids with a MW of 18.6 kDa.
TRAIL	6215 (100ug)	It is a protein that functions as a ligand to induce the process of apoptosis by binding death receptors DR4 and DR5.

Caspase Target	Fluorescent Label	Caspase Detection Kit	Catalog Numbers
Poly Caspases; apoptosis	FAM (green)	FAM-FLICA Poly Caspase assay kit	91 (25 tests); 92 (100 tests)
	SR (red)	SR-FLICA Poly Caspase assay kit	916 (25 tests); 917 (100 tests)
	660 Dye (far-red)	FLICA 660 Poly Caspase assay kit	9120 (25-50 tests)
Caspase 1 (ICE); pyroptosis, inflammation	FAM (green)	FAM-FLICA Caspase 1 assay kit	97 (25 tests); 98 (100 tests)
	660 Dye (far-red)	FLICA 660 Caspase 1 assay kit	9122 (25-50 tests)
Caspase 2; regulator	FAM (green)	FAM-FLICA Caspase 2 assay kit	918 (25 tests); 919 (100 tests)
Caspase 3/7; effectors, apoptosis (DEVD-ases)	FAM (green)	FAM-FLICA Caspase 3/7 assay kit	93 (25 tests); 94 (100 tests)
Caspase 6; effector, apoptosis	FAM (green)	FAM-FLICA Caspase 6 assay kit	95 (25 tests); 96 (100 tests)
Caspase 8; initiator, extrinsic apoptosis pathway	FAM (green)	FAM-FLICA Caspase 8 assay kit	99 (25 tests); 910 (100 tests)
Caspase 9; initiator, intrinsic apoptosis pathway	FAM (green) SR (red)	FAM-FLICA Caspase 9 assay kit	912 (25 tests); 913 (100 tests)
		SR-FLICA Caspase 9 assay kit	960 (25 tests); 961 (100 tests)
Caspase 10; initiator, extrinsic apoptosis pathway	FAM (green)	FAM-FLICA Caspase 10 assay kit	922 (25 tests); 923 (100 tests)
Caspase 13; effector, apoptosis	FAM (green)	FAM-FLICA Caspase 13 assay kit	929 (25 tests); 930 (100 tests)

Apoptosis in primary rat hippocampal neurons

Mouse EL4 thymoma cells were treated without (A) or with (B) 100ug mitomycin for 18 hours. Blue line: unstained control; Red line: SR-DEVD-FMK inhibitor probe to detect active caspases 3&7.

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FAM-FLICA green caspase detection系列 : reagent, fixative, wash buffer, Hoechst and PI
 SR-FLICA red caspase detection系列 : reagent, fixative, wash buffer, and Hoechst
 FLICA 660 far-red caspase detection系列 : reagent, fixative, wash buffer