

Plan A

NucleoSpin RNA

BDM2017014

Product	Preps	Specification	REF
NucleoSpin RNA	10	NucleoSpin RNA Columns with Collection Tubes NucleoSpin Filters Collection Tubes (2 mL) Collection Tubes (1.5 mL) Buffers <u>RNase-free rDNase</u>	740955.10
	50	as above	740955.50
	250	as above	740955.250

試用回購，特惠中

(p.s. 10 次試用品，半價試用，可抵首購金額)


Product at a glance

- Typical yield: 14μg (10⁶ HeLa; 70μg (10⁹ bacterial cells)
- A₂₆₀/A₂₈₀: 1.9—2.1
- Typical RIN: > 9
- Fragment size: > 200 nt
- Elution volume: 40—120 μL

Sample type

- Cultured cell
- Bacteria
- Yeast
- Tissue

Protocol



優惠期限: 20171001-20171231

服務專線：0800-211667

Email : service@biogenesis.com.tw

台北 TEL：02-8792-9619

台北 FAX：02-8792-9620

Website : www.bio-genesis.com.tw

Plan B

NucleoZol



740404.200



NucleoZOL	Competitor Zol
 Sample homogenization	 Sample homogenization
 Addition of <u>water</u> + Non-toxic	 Addition of <u>chloroform</u> - Toxic
 Centrifugation at RT Precipitation of contaminants	 Centrifugation at 4°C Liquid-liquid phase separation
 Aspiration of whole supernatant with RNA + <u>Easy sampling of RNA</u>	 Aspiration of aqueous phase with RNA - <u>Inconvenient pipetting for phase separation</u> - <u>Risk of carry-over of interphase/ polar phase</u> - <u>Refrigerated centrifuge necessary</u>
 Centrifugation at RT Precipitation of RNA	 Centrifugation at 4°C Precipitation of RNA
 Washing of RNA RNA in pellet	 Washing of RNA RNA in pellet
 Resuspension of RNA + <u>No drying of RNA necessary</u> + <u>Quick and easy</u>	 Drying and resuspension of RNA - <u>Time-consuming drying of RNA</u>

Related product

Advantage of NucleoZol

Product	Preps	Specification	REF
NucleoSpin RNA Set for NucleoZol	50	NucleoSpin RNA Columns Collection Tubes Buffers	740406.50

Genomic DNA Kit 特惠中

Product	Preps	Specification	REF
NucleoSpin® Tissue	50	NucleoSpin® Tissue Columns, Collection Tubes (2 mL), buffers, Proteinase K	740952.50
NucleoSpin® DNA Stool	50	NucleoSpin® DNA Stool Columns, NucleoSpin® Bead Tubes Type A, NucleoSpin® Inhibitor Removal Columns, Collection Tubes (2 mL), Collection Tubes (2 mL, lid), buffers	740472.50

服務專線：0800-211667

Email : service@biogenesis.com.tw

台北 TEL：02-8792-9619

台北 FAX：02-8792-9620

Website : www.bio-genesis.com.tw