



MACHEREY-NAGEL


Filtration · Rapid Tests · Water Analysis · Chromatography · Bioanalysis
Filtration · Schnellteste · Wasseranalytik · Chromatographie · Bioanalytik

NucleoSpin -XS 系列

一顆細胞就能抽

微量檢體

珍貴檢體



Function	Product	Cat. No.
Clean up	NucleoSpin Gel and PCR Clean-up XS	740611.50
	NucleoSpin gDNA Clean-up XS	740904.50
	NucleoSpin RNA Clean-up XS	740903.50
DNA	NucleoSpin Tissue XS;	740901.50
	NucleoSpin cfDNA XS	740900.50
	NucleoSpin DNA FFPE XS	740980.50
RNA	NucleoSpin RNA XS;	740902.50
	NucleoSpin RNA Plus XS	740990.50
	NucleoSpin totalRNA FFPE XS	740969.50

NucleoSpin Gel and PCR Clean-up

四大特色：

- ※ 同一組 kit 可用於 PCR clen 與 gel extraction
- ※ 高回收率(>80%)，回收 DNA 片段。50 bp—20 kbp
- ※ 回收體績小 (15 uL)
- ※ 含 pH 指示測，用顏色直接判斷最佳反應條件



直覺式操作：用顏色直接判斷最佳反應條件，省時又方便！



最佳反應條件，
可得最佳結果！




此反應條件，
結果稍受影響



此反應條件，
結果受到影響



50bp 片段的回收率大於 80%



Cat. No.	產品	Rxn
740609.50	NucleoSpin Gel and PCR clean up	50

服務專線：0800-211667

台北 TEL：02-8792-9619

Email：service@biogenesis.com.tw

台北 FAX：02-8792-9620

Website：www.biogenesis.com.tw