

1:10,000,000



Cloud-Clone Corp's Loading Control Antibody

可稀釋 1 千萬倍的 internal control antibody



MW	Cytoplasm	Nucleus	Cytoskeleton
>120 kDa			
70-80 kDa		Lamin A/C (63 kDa 73 kDa)	
60-70 kDa		Lamin B1 (66 kDa) ★	
50-60 kDa			
40-50 kDa			Beta Tubulin (50 kDa) ★ Beta Actin (42 kDa) ★
30-40 kDa	GAPDH (36 kDa) ★		
20-30 kDa		PCNA (29 kDa)	
<20 kDa			




Fig 1. WB - Anti-GAPDH antibody
 Line 1: Rat lung lysate;
 Line 2: Rat heart lysate
 Line 3: Mouse brain lysate;
 Line 4: Mouse kidney lysate
 Line 5: Porcine brain lysate;




Fig 2. WB - Anti-GAPDH antibody
 Line 1: DU145 cell lysate ;
 Line 2: OS-RC-2 cell lysate;
 Line 3: T-47D cell lysate
 Line 4: HEC-1B cell lysate ;




Fig 3. WB - Anti-Beta-Actin antibody
 Line 1: DU145 cell lysate ;
 Line 2: THP1 cell lysate
 Line 3: OS-RC-2 cell lysate;
 Line 4: MDA-MB -435S cell lysate ;

Catalog	Name
CAA280Hu22	Anti-Calnexin Monoclonal Antibody
CAA280Hu01	Anti-Calnexin Polyclonal Antibody
CAA285Hu22	Anti-Histone H3 Monoclonal Antibody
CAA285Hu01	Anti-Histone H3 Polyclonal Antibody
CAA285Mi01	Anti-Histone H3 Polyclonal Antibody
CAA591Hu22	Anti-PCNA Monoclonal Antibody
CAA591Hu01	Anti-PCNA Polyclonal Antibody
CAA930Hu22	Anti-ERK2 Monoclonal Antibody
CAA930Hu01	Anti-ERK2 Polyclonal Antibody
CAB021Hu22	Anti- β -catenin Monoclonal Antibody
CAB021Hu01	Anti- β -catenin Polyclonal Antibody
CAB021Mi22	Anti- β -catenin Monoclonal Antibody
CAB021Mi01	Anti- β -catenin Polyclonal Antibody
CAB340Hu22	Anti- β Actin Monoclonal Antibody
CAB340Hu01	Anti- β Actin Polyclonal Antibody
CAB340Mi22	Anti- β Actin Monoclonal Antibody
CAB340Mi01	Anti- β Actin Polyclonal Antibody
CAB341Hu22	Anti-Actin Alpha 1, Cardiac Muscle (ACTC1) Monoclonal Antibody
CAB341Hu01	Anti-Actin Alpha 1, Cardiac Muscle (ACTC1) Polyclonal Antibody
CAB342Hu22	Anti-Actin Alpha 2, Smooth Muscle (ACTa2) Monoclonal Antibody
CAB342Hu01	Anti-Actin Alpha 2, Smooth Muscle (ACTa2) Polyclonal Antibody
CAB342Mi01	Anti-Actin Alpha 2, Smooth Muscle (ACTa2) Polyclonal Antibody
CAB357Hu22	Anti-ERK1 Monoclonal Antibody
CAB357Hu01	Anti-ERK1 Polyclonal Antibody
CAB357Mi22	Anti-ERK1 Monoclonal Antibody
CAB357Mi01	Anti-ERK1 Polyclonal Antibody
CAB839Hu22	Anti-Vinculin Monoclonal Antibody
CAB839Hu01	Anti-Vinculin Polyclonal Antibody
CAB870Hu22	Anti-Tubulin β Monoclonal Antibody
CAB870Hu01	Anti-Tubulin β Polyclonal Antibody
CAB870Mi22	Anti-Tubulin β Monoclonal Antibody
CAB870Mi01	Anti-Tubulin β Polyclonal Antibody
CAB932Hu22	Anti-GAPDH Monoclonal Antibody
CAB932Hu01	Anti-GAPDH Polyclonal Antibody
CAB932Mi01	Anti-GAPDH Polyclonal Antibody
CAB932Mi22	Anti-GAPDH Monoclonal Antibody
CAC060Hu22	Anti-PDGFRa Monoclonal Antibody
CAC060Hu01	Anti-PDGFRa Polyclonal Antibody
CAF548Hu22	Anti-Lamin B1 Monoclonal Antibody
CAF548Hu01	Anti-Lamin B1 Polyclonal Antibody
CAF550Hu22	Anti-Lamin A/C Monoclonal Antibody
CAF550Hu01	Anti-Lamin A/C Polyclonal Antibody
CAG369Hu22	Anti-ATP1a1 Monoclonal Antibody
CAG369Hu01	Anti-ATP1a Polyclonal Antibody
CAG901Hu22	Anti-VDAC1 Monoclonal Antibody
CAG901Hu01	Anti-VDAC1 Polyclonal Antibody